Mölnlycke®

Mepitel® One

WITH SAFETAC® TECHNOLOGY

Single use
Usage unique

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Soft silicone wound contact layer

Product description

Mepitel® One is a transparent, perforated and non-absorbent dressing. The perforated structure of Mepitel One allows exudate to pass into an outer absorbent dressing. This can reduce frequent dressing changes closest to the wound and allows for changing the outer dressing with minimised pain and trauma. Mepitel One can also be used as a stand-alone dressing for protection of damaged skin. The dressing has a Safetac® wound contact layer which is a unique adhesive technology. It minimises pain to patients and trauma to wounds and the surrounding skin at dressing removal.

Mepitel One consists of:

- a soft silicone wound contact layer (Safetac) on one side
- a transparent, flexible, thin and perforated polyurethane film

Intended use

Mepitel One is a wound contact layer designed for the management of a wide range of exuding wounds such as; skin tears, skin abrasions, surgical incisions, partial thickness burns, traumatic wounds, partial and full thickness grafts, radiated skin, leg and foot ulcers.

It can also be used as a protective layer on non-exuding wounds, blisters and on areas with fragile skin.

Warnings

- When using Mepitel One with NPWT, always document the number of pieces of Mepitel One used in the patient's notes. This is important to ensure that no Mepitel One is left in the wound when the dressing is changed.

Precautions

- If you see signs of infection e.g. fever or the wound or surrounding skin becoming red, warm or swollen, consult a health care professional for appropriate treatment.
- Always consult a health care professional before using Mepitel One on Epidermolysis Bullosa patients.

- When using Mepitel One on partial thickness burns with high risk of rapid granulation or after facial resurfacing, avoid placing pressure on the dressing. Lift and reposition the dressing at least every second day.
- When using Mepitel One for the fixation of skin grafts and protection of blisters, the dressing should not be changed before the fifth day post application.
- Do not use Mepitel One on patients with known sensitivity to silicone or polyurethane.
- Do not reuse. If reused performance of the product may deteriorate, cross contamination may occur.
- Sterile. Do not use if inner package is damaged or opened prior to use. Do not re-sterilise.

Instructions for use

Mepitel One can be used by health care professionals or lay-persons

Mepitel One with this symbol can help you apply and remove the dressing on skin tears. For these wounds, there is a risk of re-opening the flap when removing the dressing. Apply the dressing with the arrow pointing in the same direction as the flap. Start application from the base of the flap.

When removing the dressing, begin removal following the printed symbol, in the direction the arrow is pointing.

1. Cleanse the wound in accordance with clinical practice and dry the surrounding skin thoroughly.
2. Choose a size of Mepitel One that covers the wound and the surrounding skin by at least 2 cm. If needed, the dressing can be cut. For larger wounds, more overlap is required.
3. Remove the protective plastic layers and apply Mepitel One with the sticky side to the wound. The dressing is applied in a correct way when you can read the text printed on the dressing. Smooth the dressing in place onto the surrounding skin to get a good seal. If more than one piece of Mepitel One is used, overlap the dressings. Make sure that the holes are not blocked.
4. Apply an outer absorbent dressing, on top of Mepitel One and fixate.

Minimise overlap and wrinkles, especially when Mepitel One is applied circular or wrapped. This will reduce the risk of stasis of the wound area and limb.

For application guidance see illustrations in this leaflet.

Mepitel One can be used under compression products. Mepitel One can be used together with Negative Pressure Wound Therapy (NPWT). Mepitel One can be used in combination with gels.

When the outer absorbent dressing is saturated it should be changed with Mepitel One left in place. Mepitel One is a single use dressing.

Disposal should be handled according to local environmental procedures.

Mepitel One may be left in place for up to 14 days depending on the condition of the wound and surrounding skin, or as indicated by clinical practice. To prevent maceration, exudate should pass freely through the dressing and the holes should not be blocked. When used with NPWT, follow the clinician's recommendation regarding frequency of change.

If any serious incident has occurred in relation to Mepitel One, it should be reported to Mölnlycke Health Care.

Mepitel and Safetac are registered trademarks of Mölnlycke Health Care AB.

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CANADA Mölnlycke Health Care Inc.,
2010 Winston Park Drive, Suite 100,
Oakville, Ontario, L6H 5R7
Tel: +1 905 829 1502

CHINA Mölnlycke Health Care
(Shanghai) Co., Ltd.
RM 629, No.8 Hua Jing Road,
Wai Gao Qiao FTZ Shanghai, China
Tel: +86 10 5128 8571

CZECH REPUBLIC Mölnlycke Health Care s.r.o.
Hajkova 2747/22, 130 00 Prague 3
Tel: +420 221 890 517 (Reception)
Tel: +420 221 890 511 (Customer Service)

DENMARK Mölnlycke Health Care ApS,
Gydevang 33, DK-3450 Allerød
Tel: +45 48 16 82 68
Customer Service Tel: +80 886 810

ESTONIA Mölnlycke Health Care OÜ,
Vabaõhumuseumitee 4-67, EE-13522 Tallin
Tel: +372 671 1520

FINLAND Mölnlycke Health Care Oy,
Pitäjänmäentie 14, 3 krs, FI-00380 Helsinki
Tel: +358 201 622 300
Customer Service Tel: 0800 113 442

FRANCE Mölnlycke Health Care SAS,
13 Allée du Château Blanc
59290 WASQUEHAL
Tel: +33 320 122 555
Customer Service Tel: 0800 910 292

GERMANY Mölnlycke Health Care GmbH
Grafenberger Allee 297
DE-40237 Düsseldorf
Tel: +49 211 92 08 80
Customer Service Tel: 0800 186 21 80

HUNGARY Mölnlycke Health Care Kft
HU-1134 Budapest, Dévai u. 26-28. III. em.
Tel: +36 1 477 3080

INDIA Mölnlycke Health Care India Pvt Ltd
Universal Business Centre
Second Floor, Universal Trade Tower
Sector 49, Gurgaon Sohna Road
Gurgaon 122 018, Haryana
Tel: +91 124 4696 008

ITALY Mölnlycke Health Care s.r.l.,
Via Marsala 40/C, IT-21013 Gallarate (VA)
Tel: +39 0331 714 011

JAPAN Mölnlycke Health Care KK,
4F Concieria Nishi-Shinjuku Tower's West
6-20-7 Nishi-Shinjuku,
Shinjuku-Ku, Tokyo 160-0023
Tel: +81 3 6914 5004

KOREA Mölnlycke Health Care Korea Co., Ltd
2nd floor, Hansol Richville.
46 Wiryeseong-daero,
Songpa-gu, Seoul, 05627
South Korea
Tel: +82 2 3789 1402

LATVIA Mölnlycke Health Care AB,
Dzelzavas 120 M, Riga LV-1021
Tel: +372 6711 522

LITHUANIA Mölnlycke Health Care AB,
J. Galvydzio g. 3, LT 08236 Vilnius
Tel: +370 5 2274 58 37

NETHERLANDS Mölnlycke Health Care B.V.
Postbus 3196, 4800 DD Breda
Tel: +31 76 521 9663

NORWAY Mölnlycke Health Care AS
Postboks 6229 Etterstad, NO-0603 Oslo
Tel: +47 22 70 63 70
Customer Service Tel: 800 161 78

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Mölnlycke Health Care Ltd.
USA Mölnlycke Health Care US, LLC