

**PRODUCT DATA SHEET
GAUZE- NON-STERILE ASSORTMENT**

Product name	Gauze					
Product description	Gauze swabs and sponges – non-sterile					
Intended use	Wound cleansing, skin cleansing e.g. before a surgical operation, and absorption of exudates.					
Product Codes	Art no	Description	Size cm	Size	Pack count	Comments
Non-sterile	153000-15	Swabs, 12 ply	5x5		120/6120	
	153100-15	Swabs, 8 ply	5x5		150/7650	
	153300-15	Swabs, 8 ply	7.5x7.5		150/3600	
	153700-15	Swabs, 12 ply	7.5x7.5		120/2880	
	152400-15	Swabs, 8 ply	7x10		150/3150	
	153600-15	Swabs, 8 ply	10x10		120/2400	
	158400-15	Swabs, 12 ply	10x10		100/2000	
	158500-15	Swabs, 12 ply	10x20		50/1000	
	151600-16	Sponges		no 2	500/9000	Flat
	151700-16	Sponges		no 3	500/10000	Flat
	150700-16	Sponges		no 2	250/6500	Round
	150800-17	Sponges		no 3	250/3000	Round
	151300-16	Sponges		no 4	100/1500	Round
	151400-16	Sponges		no 5	100/1500	Round
Manufacturer	Mölnlycke Health Care AB Box 13080 SE-40252 Göteborg Sweden					
Country of origin	China					
Product classification according to MDD	Class I, Non-sterile					
CE-marked	Yes					
Sterilization method	Non-sterile					
Shelf life	NA					
Product material	Gauze	100% cotton				

This document has been printed by the PRIME system. The validity of this document cannot be guaranteed.

© Mölnlycke Health Care AB. This document is the property of Mölnlycke Health Care and must not be reproduced, disclosed to any third party or used in any unauthorised manner without written consent.

	Elastic ring for sponges	Silicone		
Environmental classification	ISO 14001	Green		
Other info of special interest	Bleached without chlorine			
Physical characteristics	Gauze quality assurance is performed in accordance to EN 14079.			
Packaging	The purpose of the dispenser box and transport box is to contain and protect material from dust, moisture and damage during transport, handling and storage.			
Storage Condition	N/A			
GS1 Codes	Art No	CON	RET	TRP
Non-Sterile	153000-15		7332430418364	7332430416100
	153100-15		7332430418371	7332430416124
	153300-15		7332430418388	7332430416148
	153700-15		7332430418401	7332430416186
	152400-15		7332430418357	7332430416087
	153600-15		7332430418395	7332430416162
	158400-15		7332430418418	7332430416209
	158500-15		7332430418425	7332430416223
	151600-16		7332551208448	7332551208431
	151700-16		7332551208615	7332551208608
	150700-16		7332551208691	7332551208684
	150800-17		7332551208943	7332551208936
	151300-16		7332551208974	7332551208967
	151400-16		7332551209001	7332551208998