

Monopolar instruments

Truly customised laparoscopic solutions



A comprehensive range of premium-quality single-use 5mm monopolar instruments with two different handle designs, providing surgeons with the ultimate choice and convenience. Available both as sterile single-packed product and for inclusion in Mölnlycke® procedure trays.

Reusable-style handle

- Ergonomic reusable-style handle ensures surgeon comfort and control
- Specific handle design provides the opportunity to convert from reusable instruments to single-use equivalents without compromise
- Simple intuitive ratchet mechanism
- Ratchet mechanism can be 'turned off' – so you can use the same instrument with or without ratchet – for maximum convenience and flexibility
- Robust, wide-opening, dual-articulating jaw tips provide optimal application during a procedure
- Extensive range of jaw-tip configurations providing complete surgical offering
- Fully-insulated shaft ensures effective current containment

Single-use-style handle

- Ergonomic single-use-style handle ensures surgeon comfort and control
- Simple intuitive ratchet mechanism
- Robust, wide-opening, dual-articulating jaw tips
- Extensive range of jaw-tip configurations provide complete surgical offering
- Fully-insulated shaft ensures effective current containment



Straight Grasping Forceps




Maryland Dissecting



Clinch Grasping Forceps



Babcock Grasping Forceps



Rat Tooth Grasping Forceps



Duckbill Grasping Forceps



Straight Grasping Forceps – Fenestrated



Intestinal Grasping Forceps

Single-use-style



Single-use-style handle



Single-use-style handle Maryland

Ref. No.	Description	Size	Length	Pcs/disp.
899200	Straight grasping forceps with ratchet	5mm	33cm	10
899201	Maryland dissecting forceps	5mm	33cm	10
899202	Clinch grasping forceps with ratchet	5mm	33cm	10
899203	Babcock grasping forceps with ratchet	5mm	33cm	10
899204	Rat tooth grasping forceps with ratchet	5mm	33cm	10
899205	Duckbill grasping forceps with ratchet	5mm	33cm	10
899206	Intestinal grasping forceps with ratchet	5mm	33cm	10
899207	Straight grasping forceps with ratchet – fenestrated	5mm	33cm	10

Reusable-style



Reusable-style handle



Reusable-style handle Maryland

Ref. No.	Description	Size	Length	Pcs/disp.
899208	Straight grasping forceps with ratchet control	5mm	33cm	10
899209	Maryland dissecting forceps	5mm	33cm	10
899210	Clinch grasping forceps with ratchet control	5mm	33cm	10
899211	Babcock grasping forceps with ratchet control	5mm	33cm	10
899212	Rat tooth grasping forceps with ratchet control	5mm	33cm	10
899213	Duckbill grasping forceps with ratchet control	5mm	33cm	10
899214	Intestinal grasping forceps with ratchet control	5mm	33cm	10
899215	Straight grasping forceps with ratchet control – fenestrated	5mm	33cm	10

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, P.O. Box 13080, Gamlestadvägen 3 C, SE-402 52 Göteborg, Sweden. Phone + 46 31 722 30 00
 The Mölnlycke trademark, name and logo is registered globally to one or more of the Mölnlycke Health Care Group of Companies.
 © 2017 Mölnlycke Health Care AB. All rights reserved. HQSU0323